

**Watson-Marlow Pumps Group**  
 Telephone +44 (0) 1326 370370  
 Fax +44 (0) 1326 376009  
 E-mail: [info@watson-marlow.co.uk](mailto:info@watson-marlow.co.uk)  
[www.watson-marlow.com](http://www.watson-marlow.com)




The information contained in this document is believed to be correct, but the Watson-Marlow Pumps Group accepts no liability for any errors it contains, and reserves the right to alter specifications without notice.


**WARNING**  
 This product is not designed for use in, and should not be used for, patient-connected applications.

PB0294/10

	101F - 50Hz - 4rpm, 20rpm 101F - 60Hz - 4.8rpm, 24rpm 101U - 0.06 - 2rpm 101U - 1.0 - 32rpm		+5C, +41F +40C, +104F -40C, -40F +70C, +158F
	<70dBA@1m		101F - 1.7kg (4lb) 101U - 2.2kg (4.8lb)



*David Cole*  
 David Cole, Managing director



**101F/R, 101U/R**