






A

040.3*1V.M*U = 240V/1/50/60Hz
 040.3*1V.M*E = 230V/1/50/60Hz
 040.3*1V.M*A = 110V/1/60Hz
 040.5*1V.M*U = 240V/1/50/60Hz
 040.5*1V.M*E = 230V/1/50/60Hz
 040.5*1V.M*A = 110V/1/60Hz



403U/VM2 403U/VM3 403U/VM4



Watson-Marlow Pumps Group
 Telephone +44 (0) 1326 370370
 Fax +44 (0) 1326 376009
 E-mail: info@watson-marlow.co.uk
www.watson-marlow.com



The information contained in this document is believed to be correct, but the Watson-Marlow Pumps Group accepts no liability for any errors it contains, and reserves the right to alter specifications without notice.



WARNING
 This product is not designed for use in, and should not be used for, patient-connected applications.

PB0313/10

	Max 10rpm 50rpm		+5C, +41F
	<70dBA@1m		+40C, +104F
	403U/VM2 - 2.5kg 403U/VM3 - 2.5kg 403U/VM4 - 2.5kg		-40C, -40F
			+70C, +158F



403U/VM2, 403U/VM3, 403U/VM4



David Cole
 David Cole, Managing director

2006/42/EC, 2004/108/EC,
 2006/95/EC, BS EN60529 (IP21),
 2006/42/EC

